



Full Face Mask (F6)

To redefine the comfort of full face therapy
family cares



Full Face Mask (F6)

Main features




Under-the-nose Cushion

The innovative under-the-nose cushion brings liberty for the nasal bridge and meanwhile minimizes the claustrophobia.

Circular Dispersion Vent

The circular dispersion vent keep the noise as low as 17dB (A) and meanwhile protect your bed partner from being disturbed.



Modular Main Frame

The main frame is compatible with all cushion sizes.


Crown-style Adjustable Headgear

The adjustable headgear is fit for most of the head sizes. Wider range, more inclusiveness.



Technical Specifications

Pressure–Flow Curve



Pressure (cmH ₂ O)	4	10	12	17	24	30
Flow Rate (L/min)	22	37	41	49	59	67

Dead Space Information

Dead space is the empty volume of the mask up to the swivel. The dead space of the mask varies according to cushion sizes but is less than 180ml.

Therapy Pressure

4 to 30 cmH₂O

Resistance

Drop in Pressure measured (average for 3 sizes)
at 50 L/min: 0.15 cmH₂O, at 100 L/min: 0.5 cmH₂O

Inspiratory and Expiratory Resistance

The inspiratory resistance of the mask (in combination with the Non-Rebreathing Valve) is 1.8 cmH₂O at 50 L/min.
The expiratory resistance of the mask (in combination with the Non-Rebreathing Valve) is 2.0 cmH₂O at 50 L/min.

Sound

DECLARED DUAL-NUMBER NOISE EMISSION VALUES in accordance with ISO 4871. The A-weighted sound power level of the mask is 25 dBA, with uncertainty 3 dBA. The A-weighted sound pressure level of the mask at a distance of 1 m is 17 dBA, with uncertainty 3 dBA.



F5 Full Face Mask

Simplify the Full Face Therapy

F5 do renovates your impression of full face mask with slim frame, groove- fit cushion and easy clip. All of the components above enable the patient to enjoy the ease of full face therapy.



Slim Frame

Transparent and slim main frame minimize the effect to your vision.

Groove-fit Cushion

Single layer of soft cushion with intelligent groove design fits different face shapes.

Easy Clip

Classic clip helps to lock or release as easily as possible.

Ingenious Swivel

Detachable swivel design separates mask from tubing conveniently when needed.

Product Code

F5 Full Face Mask

10030268 (L)

10030269 (M)

10030270 (S)



Frame Assembly

60380294



Cushion Assembly

60380295 (L) 60380296 (M) 60380297 (S)



BMC Medical Co., Ltd.

Add: Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA

Tel: +86-10-51663880

Fax: +86-10-51663880 ext.810

en.bmc-medical.com



F5A Full Face Mask

Enlarge your vision by efficient face mask

F5A is one of the newest full face series. It is capable of auto adaptation to the nose bridge. Simple structure and light weight provide you with more comfortable and effective therapy.



Auto-adaptation

The groove on the cushion is able to adapt itself to the nose bridge, maintaining low air leakage and limited face pressure.

Liberty of vision

Based on the forehead frame-free design, users can embrace more vision with ease.

Anti-asphyxia valve

Under the premise of security, the position of anti-asphyxia valve is adjusted in order to achieve more functions on a smaller elbow.

More comfort

Due to the single-layer transitional design, the cushion touches the nose and face softly.

Product Code

F5A Full Face Mask

10030322 (L)

10030323 (M)

10030324 (S)



Replacement Headgear for F5A



Replacement Cushion for F5A

- Large size
- Medium size
- Small size





BMC Medical Co., Ltd.

Room 110, Tower A Fengyu Building, No.115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA

Tel: +86-10-5166 3880 / Fax: +86-10-5166 3880 ext. 810 / en.bmc-medical.com

Not applicable to some markets, such as China, the U.S., the U.K., Germany and Japan

ST-SC-02-21-12-002-V1.0